Code No. 292-81101 (For 50 tests)

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FUJIFILM Wako Chemicals (Hong Kong) Limited

Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong
Telephone : + 852-2799-9019
wkhk.info@fujifilm.com

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<https://labchem-wako.fujifilm.com>

Code No. 292-81101 (For 50 tests)

DNA Extractor™ Kit for Residual DNA, CP Method (Sodium Iodide Method)

This product is in compliance with sodium Iodide method and designed for use in extracting residual host cell DNA in biological samples such as vaccines and biopharmaceuticals.

It offers high recovery of DNA even from samples containing only a very small amount of DNA. The entire extraction step is completed with a single tube for 60-90 minutes. The extracted DNA can be quantified by qPCR.

The standard protocol provided in this kit is suitable for extracting DNA from raw material samples (eg. cell lysates, or solutions containing high concentrations of proteins or lipids). At the same time, we also provide simplified steps for extracting DNA from purified samples (eg. purified proteins and serum). The differences between the two operating methods are shown in the table below.

	Sample of interest	Pretreatment step of sample	Number of washing steps	"Washing Solution A, CP" in the kit
Standard Protocol	Crude	Required	2 times	Used
Simplified Protocol	Clean	Not required	1 time	Not used

[Features]

- In compliance with sodium Iodide method
- High quality and recovery
- Completed in a single tube for 60-90 minutes.

[Principle of sodium Iodide method for DNA extraction]

An chaotropic reagent, Sodium Iodide, and an anionic detergent participate in solubilization of proteins and lipids contained in biological samples. After addition of 2-propanol to the mixture, nucleic acids are co-precipitated with glycogen as a carrier, while other components remain soluble in the solution phase.

[Kit contents]

Code No.	Components	Size
299-81111	Sodium Iodide Solution, CP	26 mL × 1
296-81121	Sodium N-Lauryl Sarcosinate Solution, CP	1.2 mL × 1
293-81131	Washing Solution A, CP	42 mL × 1
290-81141	Washing Solution B, CP	40 mL × 2
297-81151	Glycogen Solution, CP	0.1 mL × 1

[Storage]

Store at 2-10°C in the dark.

[Additional materials required for 50 tests]

<Reagents>

- 1) 2-Propanol 20 mL
- 2) Sterile distilled water 6 mL

<Instruments>

- 1) Microcentrifuge
- 2) Vortex mixer
- 3) 50×2 mL centrifuge tube
- 4) Micropipette
- 5) Micropipette tips
- 6) Heating block

—5/8—

[Precautions for Use]

- The color of "Washing Solution A, CP" may be slightly different between product lots, e.g. faint yellow or colorless. It does not affect quality of DNA extraction.
- If precipitate is observed in "Sodium Iodide Solution, CP" or "Sodium N-Lauryl Sarcosinate Solution, CP", warm it or them at 50°C until the precipitate is no longer visible.
- Use DNA-free, DNase-free and sterilized reagents and Instruments.
- DNA concentration cannot be accurately measured by absorbance if sodium Iodide remains in the extracted DNA. Measure the amount of DNA by qPCR

[Standard Protocol]

<Preparation of Reagents>

1. Solution I

Add 65 μL of "Glycogen Solution, CP" into the bottle of 26 mL of "Sodium Iodide Solution, CP". Mix it with a vortex mixer.

Note : Solution I can be stored at 2-10°C in the dark for approximately 5 months.

2. Solution II

Add 2 μL of "Glycogen Solution, CP" into the bottle of 40 mL of "Washing Solution B, CP". Mix it with a vortex mixer immediately.

Notes :

- It does not affect quality of DNA extraction if white precipitates appear in Solution II.
- Solution II can be stored at 2-10°C in the dark for approximately a week at least. If preparing a small volume of Solution II, use required amount of "Glycogen Solution, CP" and "Washing Solution B, CP". For example, add 1 μL of "Glycogen Solution, CP" into 20 mL of "Washing Solution B, CP" in a sterilized tube.

<Preparation of Samples>

■ In case that protein concentration in the sample solution is less than 2 mg/mL.

1. Add the reagents below into the sample solution to the final concentration.

Reagents	Final concentration
Sodium Dodecyl Sulfate	0.10%
Dithiothreitol	50 mmol/L

2. Warm it at 55°C for 30 minutes in a heating block.

■ In case that protein concentration in the sample solution is more than 2 mg/mL.

1. Add the reagents below into the sample solution to the final concentration.

Reagents	Final concentration
Sodium Dodecyl Sulfate	0.10%
Dithiothreitol	50 mmol/L
Sodium Chloride	150 mmol/L
EDTA	1 mmol/L

2. Add 20 μg of Proteinase K per 1 mg of protein in the sample solution in to the sample solution.

3. Adjust pH to be approximately 7.5 with Tris.

4. Warm it at 55°C for an hour at least in a heating block.

<DNA Extraction Procedure>

1. Dispense 400-500 μL of sample solution to a 2 mL blank centrifuge tube.
2. Add 20 μL of "Sodium N-Lauryl Sarcosinate Solution, CP" into the tube and mix it with a vortex mixer.
3. Add 500 μL of Solution I into the tube and mix it with a vortex mixer.
4. Warm the tube at 40°C for 15 minutes in a heating block.
5. Remove the tube from the heating block. Add 900 μL of 2-Propanol into the tube and mix it with a vortex mixer.
6. Leave the tube at room temperature for 15 minutes.
7. Centrifuge the tube with 10,000 × g at room temperature for 15 minutes. A faint white pellet appears in the tube.
8. Remove supernatant from the tube by using a pipette or decanting carefully until the liquid volume in the tube is less than approximately 100 μL.

—6/8—

Notes :

- Be careful not to remove the pellet from the tube.
- Centrifuge the tube with $10,000 \times g$ at room temperature for 15 minutes again if the pellet floats.
- In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.

9. Add $800 \mu\text{L}$ of "Washing Solution A, CP" into the tube and mix it by inversion mixing. Ensure that the pellet is detached from the tube wall.
10. Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes.
11. Remove supernatant from the tube by using a pipette or decanting carefully until the liquid volume in the tube is less than approximately $100 \mu\text{L}$.

Notes :

- Be careful not to remove the pellet from the tube.
- Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes again if the pellet floats.
- In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.

12. Add 1.5 mL of Solution II into the tube and mix it by inversion mixing.
13. Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes.
14. Remove as much supernatant as possible from the tube by using a pipette or decanting carefully.

Notes :

- Be careful not to remove the pellet from the tube.
- Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes again if the pellet floats.
- In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.

15. Dry the pellet by warming it at 60°C in a heating block, vacuum dryer or natural drying.

Notes :

- Too dry pellet is difficult to be dissolved. Stop drying the pellet when it still is wet.
- Much residual liquid in the tube may cause low qPCR efficiency because Solution II contains ethanol, which can inhibit qPCR.
- The pellet contains DNA and glycogen.

16. Dissolve the pellet with sterile distilled water or buffer, and proceed with qPCR.

[Simplified Protocol]

This protocol is a DNA extraction method for purified proteins and serum samples that are relatively low in impurities.

Sample pretreatment is omitted and the number of washing step is one. This reduces the operating process and operating time.

<Preparation of Reagents>

1. Solution I (Simplified version)

Add 6 mL of sterile distilled water, 1 mL of "Sodium N-Lauryl Sarcosinate Solution, CP" and $65 \mu\text{L}$ of "Glycogen Solution, CP" into the bottle of 26 mL of "Sodium Iodide Solution, CP". Mix it with a vortex mixer.

Note : Solution I (Simplified version) can be stored at $2-10^\circ\text{C}$ in the dark for approximately 5 months.

2. Solution II

Add $2 \mu\text{L}$ of "Glycogen Solution, CP" into the bottle of 40 mL of "Washing Solution B, CP". Mix it with a vortex mixer immediately.

Notes :

- It does not affect quality of DNA extraction if white precipitates appear in Solution II.
- Solution II can be stored at $2-10^\circ\text{C}$ in the dark for approximately a week at least. If preparing a small volume of Solution II, use required amount of "Glycogen Solution, CP" and "Washing Solution B, CP". For example, add $1 \mu\text{L}$ of "Glycogen Solution, CP" into 20 mL of "Washing Solution B, CP" in a sterilized tube.

<DNA Extraction Procedure>

1. Dispense $100 \mu\text{L}$ of sample solution to a 2 mL or 1.5 mL blank centrifuge tube.

2. Add $300 \mu\text{L}$ of Solution I (Simplified version) into the tube and mix it with a vortex mixer.
3. Warm the tube at 60°C for 15 minutes in a heating block.
4. Remove the tube from the heating block. Add $400 \mu\text{L}$ of 2-Propanol into the tube and mix it with a vortex mixer.
5. Leave the tube at room temperature for 15 minutes.
6. Centrifuge the tube with $10,000 \times g$ at room temperature for 15 minutes. A faint white pellet appears in the tube.
7. Remove supernatant from the tube by using a pipette or decanting carefully until the liquid volume in the tube is less than approximately $100 \mu\text{L}$.

Notes :

- Be careful not to remove the pellet from the tube.
- Centrifuge the tube with $10,000 \times g$ at room temperature for 15 minutes again if the pellet floats.
- In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.

8. Add 1 mL of Solution II into the tube and mix it by inversion mixing. Ensure that the pellet is detached from the tube wall.
9. Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes.
10. Remove as much supernatant as possible from the tube by using a pipette or decanting carefully.

Notes :

- Be careful not to remove the pellet from the tube.
- Centrifuge the tube with $10,000 \times g$ at room temperature for 5 minutes again if the pellet floats.
- In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.

11. Dry the pellet by warming it at 60°C in a heating block, vacuum dryer or natural drying.

Notes :

- Too dry pellet is difficult to be dissolved. Stop drying the pellet when it still is wet.
- Much residual liquid in the tube may cause low qPCR efficiency because Solution II contains ethanol, which can inhibit qPCR.
- The pellet contains DNA and glycogen.

12. Dissolve the pellet with sterile distilled water or buffer, and proceed with qPCR.

[FAQ]

Q. Peak wavelength appeared at 230 nm when extracted DNA concentration was measured by absorbance.

A. Sodium Iodide remains in the extracted DNA. Measure the amount of DNA by qPCR

Q. DNA recovery rate was too low.

A. A part of DNA pellet may have been removed at step 7 or 10. Remove the supernatant with a pipette carefully.

<Reference>

Ishizawa, M., Kobayashi, Y., Miyamura, T. and Matuura, S. (1991) Simple Procedure of DNA isolation from human serum, *Nucleic Acids Res.*, **19**, 5792.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964

FUJIFILM Wako Chemicals (Hong Kong) Limited

Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong
Telephone : + 852-2799-9019
wtkh.info@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation

Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road, Guangzhou, 510095, China
Telephone : + 86-20-8732-6381(Guangzhou)
Telephone : + 86-21-6288-4751(Shanghai)
Telephone : + 86-10-6413-6388(Beijing)
wkgz.info@fujifilm.com

<https://labchem-wako.fujifilm.com>